GI Genius™ software

User Manual

Version 1.1 – 12th April 2023

REF CB1708-MN03-EU: User Manual of GI Genius™ software.

For software installed at a user site located in at least one of the EU Member States, the following information does apply.

The present User Manual of the GI Genius™ software also corresponds to the electronic instructions for use (eIFU) according to European Regulation 2021/2226/EU.

This eIFU refers to the medical device software: GI Genius™ software. Basic UDI-DI: 805506011SWDCBB.

This eIFU is available in PDF format. PDF format is a widely used open file format for documents that are designed to be read-only and maintain the same appearance across multiple devices and platforms. They can be read using various free software, which can be downloaded from the Internet and installed on PCs or mobile devices.

Whenever you may require a paper copy of the present document, it will be shipped to you within seven calendar days from receipt of your request, at no additional cost. You can contact Linkverse or your authorized local distributor to request a paper copy of this eIFU.

This eIFU is identified with a REF code, an issue date, and a version. In case a newer version will be released for your software, it will be made available on Linkverse website for your download, at no additional cost.

Linkverse or your authorized local distributor will inform you if a newer version of this eIFU is available for your software.

Manufacturer contact details:

Linkverse S.r.l.

Via Ostiense 131/L, 00154 Rome (Italy)

Website: <https://linkverse.com/>

E-mail: info@linkverse.com

Tel: (+39) 06 90287198

Important Information – Please Read Before Use

If the product is used outside the region, it may not operate as specified in the specifications. Please read this Instructions for Use carefully before using the product to ensure that you use it safely, effectively and for proper maintenance. Please retain this manual for future reference.

Index

- 1 Symbols and marking..... 3
- 2 Introduction 3
 - 2.1 Safety symbols..... 3
 - 2.2 Terms and abbreviations..... 3
 - 2.3 Document purpose..... 3
 - 2.4 Software version..... 4
- 3 Description 4
 - 3.1 Intended purpose and indications of the software 4
 - 3.2 User qualification..... 5
 - 3.3 GI Genius™ software compatibility 5
 - 3.4 Performance characteristics..... 5
 - 3.5 Cybersecurity measures..... 5
- 4 Operating the software with GI Genius™ hardware 5
 - 4.1 Installation and updates 5
 - 4.2 Operating the software..... 6
 - 4.3 Turn ON/OFF the GI Genius™ device..... 6
 - 4.4 GI Genius™ software information output and user interface 7
 - 4.5 Control user interface 9
 - 4.6 GI Genius™ hardware By-Pass and Emergency Shut Down 12
- 5 Operating the software with third-party hardware..... 12
 - 5.1 Minimum requirements..... 12
 - 5.2 Cybersecurity requirements 12
 - 5.3 Safety requirements 12
 - 5.4 Software operability..... 13
- 6 Troubleshooting 13
- 7 Support..... 13
- 8 Annex 1: GI Genius™ software VP compatibility list 13

1 Symbols and marking

	Legal Manufacturer		Unique Device Identifier of the medical device
	Consult the instructions for use		The product is a medical device
	QR code with website link for eIFU		Quantity in package
	Caution		Catalogue code
	Waste for separate collection according to Directive 2012/19/EU (WEEE)		Only for professional use
	CE marking of conformity to Regulation 2017/745/EU		Fragile
	Date of manufacture		Keep dry

2 Introduction

2.1 Safety symbols

This manual and this software use the safety symbols below. The symbols denote critical information. Please read them carefully.

DANGER	Imminent hazardous situation that will result in death or serious injury, if not avoided
WARNING	Potential hazardous situation that could result in death or serious injury, if not avoided
CAUTION	Potential hazardous situation that may result in minor or moderate injury, if not avoided. Can alert about unsafe practices or potential damage
NOTE	Useful information

2.2 Terms and abbreviations

The following table summarizes terms and abbreviations used in this document.

Term	Description
FoV	Field of View
MD	Medical Device
MED	Main Endoscopy Display
VP	Endoscopy Video Processor

2.3 Document purpose

This document describes the essential information on using the GI Genius™ software in a safe and effective way. Please carefully review this manual and all the equipment manuals involved in the procedure and use the software as instructed.

WARNING	Failure to follow the instructions in this manual may result in damage to and/or malfunction of the software or other equipment involved in its use. Do not use this software for any other purpose than that for which it has been intended.
CAUTION	This document does not contain any information on how to perform endoscopy procedures, please refer to the appropriate manuals for this information.

Keep this and any related document in a safe and accessible location. If you have questions or comments about this document, please contact helpdesk service section 7.

2.4 Software version

This document relates to software release v.3.1.0 and above (e.g. v.3.1.x), if not differently specified in the text. Please be sure that the software version is such before: version can be displayed by pressing the “Confirm” button on the front panel of the GI Genius™ hardware (see below) or the “Enter” key on a third-party hardware.

 “Confirm” button on the GI Genius hardware

The interface of the software does show the software label identified below (refer to **Figure 10** for an example of the appearance of the whole window):

Figure 1: Software label shown on the lower part of the user interface window

3 Description

3.1 Intended purpose and indications of the software

The GI Genius™ software is an artificial intelligence-based medical device that has been trained to process colonoscopy images containing regions consistent with colorectal lesions like polyps, including those with flat (non-polypoid) morphology. The GI Genius™ software is intended to be used by trained clinicians as an adjunct to white-light colonoscopy for the purpose of highlighting regions suspected to have visual characteristics consistent with different types of mucosal abnormalities (e.g., colorectal polyps). The target population is represented by persons undergoing colonoscopy procedures.

If Characterization support is enabled, a polyp detected and highlighted by the GI Genius™ software is consistently framed in white-light video colonoscopy, based on the visual characteristics of the detected polyp, the GI Genius™ software provides an estimation of the possible polyp histology.

WARNING	The GI Genius™ software is intended to be used as an adjunct to colonoscopy procedures and is not intended to replace endoscopist assessment or histopathological sampling.
WARNING	The whole colonoscopy video and the regions highlighted by the GI Genius™ software must independently be assessed by the endoscopist, with all the available and obtainable information, without relying primarily on GI Genius™ software output.
WARNING	The GI Genius™ software does not prescribe any clinical management action regarding the detected and highlighted polyps. The endoscopist must take appropriate actions according to the standard clinical practice.

The GI Genius™ software does not make any elaboration or alteration of the colonoscopy video streaming, it only overlays graphical markers.

The GI Genius™ software can be used in combination with the GI Genius™ hardware or other compatible hardware whose minimum requirements are specified in section 5.1.

3.2 User qualification

Users shall have been properly trained on GI Genius™ software usage and safety.

Users shall be expert clinicians on lower gastrointestinal endoscopy procedure and the training shall be based on the contents of this user manual.

3.3 GI Genius™ software compatibility

The GI Genius™ software is intended to be used in combination with the GI Genius™ hardware.

The GI Genius™ hardware are: GI Genius™ Intelligent Endoscopy Module and GI Genius™ Modules. Please refer to the user manual of corresponding hardware for further information.

Use of GI Genius™ software with third-party hardware is possible only if such hardware does meet the minimum requirements specified in section 5.1.

Please refer to Annex 1 of this user manual for the list of Video Endoscopy system compatible with the software.

DANGER	Using incompatible equipment can result in patient injury or equipment damage and makes it impossible to obtain the expected functionality.
CAUTION	The usage of OVERSCAN display functionality could result in hiding of GI Genius™ markers.

3.4 Performance characteristics

The GI Genius™ performances resulting from an internally established benchmark representing the behavior during screening and surveillance colonoscopy is summarized below. Device comparison to endoscopists has been assessed by comparing GI Genius™ real performances to a panel of 5 endoscopists, each endoscopist with an experience of at least 500 colonoscopies.

GI Genius™ comparison to endoscopists – fraction of polyps detected earlier

- Software version 3.0.0: anticipates the average endoscopist in detecting polyps on 86.5% (95% CI: 82.9% to 90.2%) of them, by considering repeated marking of the same target as a single statistical event.

Detection False Positives

- Software version 3.0.0: 1.28% (95% CI: 1.11% to 1.47%) of the video frames showing no polyps contains false positive detections (False Positive Rate), by means of a logistic regression mixed model.

The GI Genius™ performances related to the characterization function:

Characterization Accuracy

- Software version 3.0.0: non-experts with the help of GI Genius feature an accuracy of 81.25% (95% CI: 75.73% to 85.98%), while experts alone feature an accuracy of 77.69% (95% CI: 71.91% to 82.77%). Therefore, non-experts with the help of GI Genius are noninferior to experts alone.

Negative predictive Value (or NPV)

- Software version 3.0.0: 97.6% (95% CI, 94.1%–99.1%) for adenomatous histology of diminutive rectosigmoid polyps, as resulting from an independent study (doi.org/10.1016/j.cgh.2022.04.045).

Performances of software version 3.0.0 are the same of successive updated versions, if not differently specified in the present section.

3.5 Cybersecurity measures

The GI Genius™ software does not store any data. Analysis of images is performed in real-time during the procedure and no archive nor database is present.

The GI Genius™ software does not need connections to an IT network to operate.

The GI Genius™ software does run on an encrypted mass storage.

Do not attempt to download or update the GI Genius™ software yourself. Please contact your local distributor, who will update the software for you.

The GI Genius™ software does not provide user accounts.

Store the hardware where the GI Genius™ software is installed in areas accessible to authorized personnel only. GI Genius™ software shall be used by physicians only.

4 Operating the software with GI Genius™ hardware

Please refer to user manual of the GI Genius™ hardware for any information related to the hardware.

4.1 Installation and updates

The installation and any following updates of the GI Genius™ software is operated by the authorized local distributor or directly by Linkverse. As user, you are not required to install or update the software yourself.


After installation or update of version, the software does execute controls in the initialization phase to check the correctness of the procedure, so that if a problem is encountered, the software will provide feedback on its status.

In case of installation on a GI Genius™ Intelligent Endoscopy Module, the software version already installed shall be 3.0.2 or above. Please refer to the user manual of the GI Genius™ Intelligent Endoscopy Module for any information related to such device.

For any further information or in case of doubt, please refer to the contacts in section 7.

4.2 Operating the software

When the GI Genius™ software is operated with the GI Genius™ hardware, the following functions are available:

1		<p>Stand-by/ power light indicator Press it to go into GI Genius™ hardware stand-by mode or to end the GI Genius™ hardware stand-by mode</p> <p>NOTE When the stand-by mode is off the green indicator light is solid.</p>
2		<p>Overlay activation/deactivation Press it to enable and disable GI Genius™ overlay over the main endoscopy display</p> <p>NOTE When the overaly activation/deactivation indicator light is solid, the software is operating correctly</p> <p>NOTE When the indicator light is blinking contact the support (see section 7).</p>
3	   	<p>Keyboard To allow user navigation into the menu</p> <p>Arrow keys to allow user to move up/down in the menu</p> <p>Arrow keys to allow user to move left and right in the menu</p> <p>Confirm button to allow user to open the menu and confirm actions</p> <p>Back button to allow user to close the menu</p>


Figure 2: GI Genius hardware Front Panel

4.3 Turn ON/OFF the GI Genius™ device

- Press the power supply switch (I/O) to enable (I) the power grid connection;
- The device is in stand-by mode;
- Press the stand-by button to turn on the device;
- The device emits three beeps;
- The field of view is clearly identified by four corners, and a blinking green square indicator appears on the connected MED to state that the device is ready (**Figure 3**);
- Press the stand-by button again to turn off the device;
- The device is in stand-by mode;
- Press the power supply switch (I/O) to disable (O) the power grid connection.

CAUTION	After the device has been turned off, wait at least 10 seconds before turning it back on again. A fuse may be blown if you turn the power on and off repeatedly within a short period of time. Do not use a sharp object such as the tip of a pen to press the buttons on the operation panel.
----------------	--

4.4 GI Genius™ software information output and user interface

The GI Genius™ software is a video processing software that generates as output a video on the main endoscopy display that contains the original live video together with superimposed markers that will appear when a polyp is detected and will not be visible when no detection occurs.

The aim of the markers is to highlight polyps contained in the input video stream. Markers appear as green squares and are accompanied by a short sound.

The four corners identify the Field of View and their color is representative of software status. The software status indicator is shown on the bottom-right corner, providing information about the status of the software and/or its hardware support.

	Blinking green indicator is shown when the software is activated and is fully functional (Figure 3).
	Solid red indicator is shown when a software or hardware malfunction occurs (Figure 3).
	Solid red strikethrough indicator is shown when a software or hardware malfunction occurs and the software overlay is deactivated (Figure 4).
	Solid grey strikethrough indicator is shown when the overlay is deactivated (Figure 4).
	Solid grey indicator is shown when the menu is opened and the overlay is activated (Figure 5).

IF CHARACTERIZATION SUPPORT IS ENABLED, THE BELOW TAGS ARE SHOWN:


	The purple tag “adenoma” is shown when the system predicts a possible adenoma histology (Figure 6).
	The green tag “non-adenoma” is shown when the system predicts a possible non-adenoma histology (Figure 7).
	The tag “analyzing” is shown when the system is collecting the information needed to provide a possible histology prediction (Figure 8).
	The tag “no-prediction” is shown when the system is not confident enough to provide a possible histology prediction (Figure 9).

Figure 3: Example of user interface with marker, FoV corners and green status indicator

Figure 4: Example of user interface with FoV corners, overlay deactivated and red status indicator


Figure 5: Example of user interface when the menu is opened


Figure 6: Example of user interface with marker characterization (adenoma)

Figure 7: Example of user interface with marker characterization (non-adenoma)


Figure 8: Example of user interface with marker characterization (analyzing)




Figure 9: Example of user interface with marker characterization (no-prediction)

4.5 Control user interface

The GI Genius™ software user control interface (**Figure 5**) allows the user to perform the following actions:

1. Display the status and settings;
2. Beep volume regulation;
3. Field of view settings and software overlay activation/deactivation;

1		<p>On GI Genius™ hardware keyboard press the confirm button to open the menu.</p>
<p>The GI Genius™ hardware status (Figure 10) will be displayed containing the following information:</p> <ul style="list-style-type: none"> • CPU load; • GPU load; • Hardware temperature; • RX (Input) channel status indicator; • TX (Output) channel status indicator. 		
		
<p>Figure 10: Example of Status and Settings - Hardware status</p>		
<p>NOTE if the RX channel is disconnected during the endoscopy procedure the RX indicator is grey (Figure 11)</p>		
		
<p>Figure 11: Example of Status and Settings – Hardware status with RX channel disconnected</p>		
2		<p>Press the down arrow button to move to the Volume tab (Figure 12).</p>



Figure 12: Example of System Status and Settings – Volume tab

	<p>Press the confirm button to activate/deactivate the Volume.</p>
	<p>With Volume activated, press the right arrow button to select the Volume bar (Figure 12)</p>
	<p>Press the up/down arrow keys to regulate the Volume.</p>
	<p>NOTE When the volume is disabled it is clearly identified in the GI Genius™ software menu. Please refer to Figure 13.</p>
	 <p>Figure 13: Example of Status and Settings – Volume disabled</p>
	<p>Press the left arrow key to come back to the main menu tabs (Figure 12).</p>
<p>3</p> 	<p>Press the down arrow button to move to the Field View tab (Figure 14).</p>


Figure 14: Example of System Status and Settings –Field View selected

Press the confirm button **to activate/deactivate** the software overlay.


Figure 15: Example of System Status and Settings –Field view tab with overlay deactivated

NOTE Please refer to **Figure 15** for user interface when the overlay is deactivated.

NOTE Please refer to **Figure 14** and **Figure 16** for user interface when the overlay is activated.


Figure 16: Example of System Status and Settings – Field view tab with overlay activated

Press the right arrow button to select the tab **Field View**.

Press the up/down arrows keys to change the Field View Settings.

NOTE Please refer to **Figure 16** for user interface when the Field View selected is “Small”.

	NOTE Please refer to Figure 14 for user interface when the Field View selected is “Large”.
--	--

4		Press the back button to close the menu .
----------	---	---

4.6 GI Genius™ hardware By-Pass and Emergency Shut Down

CAUTION	<p>If the hardware is in stand-by or powered off, a by-pass feature that directly connect video input to video output is automatically activated, so the video directly and only flows through the hardware to not loose the video stream to be visualized on the main endoscopy display.</p> <p>If the endoscopic image disappears, the image freezes and cannot be restored, or there are interferences in the video signal:</p> <ul style="list-style-type: none"> - Put the hardware in stand-by mode with the front switch or shut it down with the back switch (I/O) to disable (O) the power grid connection; - If the video image still cannot be restored, disconnect the hardware from the endoscopy video processor and the main endoscopy display; - Reconnect the hardware to the endoscopy video processor and the main endoscopy display; - If the video image still cannot be restored, follow the procedures described in the endoscopy unit manual.
----------------	---

5 Operating the software with third-party hardware

User shall refer to the information accompanying the third-party hardware for proper use and handling of the device.

5.1 Minimum requirements

The GI Genius™ software can be installed and operated on third-party hardware, which are independently certified and registered as applicable by the third-party organization and which shall meet the following minimum requirements:

- CPU: Intel Core i7 or AMD Ryzen 7
- RAM: 8 Gb
- GPU: NVIDIA GeForce RTX 2080 or NVIDIA Quadro RTX 3000
- Video I/O: Deltacast DELTA 3G-elp-key 11 or Deltacast DELTA 12G-elp-key 11
- OS: Ubuntu 18.04
- Mass storage: 250 Gb SSD
- Keyboard: USB

5.2 Cybersecurity requirements

- Implemented encryption of mass storage
- Restricted system user access
- Disabled networking
- Disabled user copying, downloading, and updating of GI Genius™ software
- Implemented audit trail of system events
- Implemented antitampering mechanism in the hardware.

5.3 Safety requirements

DANGER	In case of problems to the images on screen, such as loss of image continuity or disturbances, it is necessary for the third-party hardware to provide a by-pass function or emergency shutdown function to allow the user to never lose control of the images on screen.
---------------	---

The hardware shall comply to requirements of IEC 60601-1 for basic safety and essential performance of medical electrical equipment.

The hardware shall comply to requirements of IEC 60601-1-2 for electromagnetic disturbances.

The hardware shall be equipped with medical grade power supply.

The hardware shall have a shutdown mechanism or equivalent feature to allow by-pass of the video stream directly to the display, in case of fault.

The hardware shall be tested against over-temperature, resistance to heat, mechanical failure conditions and provide shutdown mechanisms or equivalent solutions to protect safe use.

The hardware shall inform about consequences of drop of liquids.

The hardware shall comply with applicable requirements on waste of electrical and electronic equipment and on restriction of the use of certain hazardous substances in electrical and electronic equipment.

The hardware shall be tested against allowed limits of emission of noise levels.

5.4 Software operability

When installed on a third-party hardware, the GI Genius™ software is operable through the following keyboard keys:

Keyboard	
To allow user navigation into the software menu	
	Spacebar button - Overlay activation/deactivation Press it to enable and disable the GI Genius™ software AI overlay on screen
	Arrow keys to allow user to move up/down in the menu
	Arrow keys to allow user to move left and right in the menu
	Enter button to allow user to open the menu and confirm actions
	Backspace button to allow user to close the menu

6 Troubleshooting

Symptom	Possible Cause	Solution	Reference
GI Genius™ software output does not appear on MED	Software failure	Contact Support	section 7
	Software Cybersecurity failure	Contact Support	section 7
The screen remains frozen	Software problem	GI Genius™ software immediate power off.	section 4.3 section 4.6 section 7
Black screen	Software problem	GI Genius™ software immediate power off.	section 4.3 section 4.6 section 7
	Software failure	Contact Support	section 7

NOTE	For any other problem please contact the help desk as provided in section 7.
-------------	--

7 Support

The helpdesk service provided to the users of the GI Genius™ software and hardware is available at the following references:

Email address

serviceandrepair.uk@medtronic.com

gi-support@medtronic.com

NOTE	In the event that a serious incident has occurred related to device use, immediately report the event to the manufacturer, competent authorities, and any other regulators as required.
-------------	---

8 Annex 1: GI Genius™ software VP compatibility list

Below is a list of VPs compatible with GI Genius™ software, through connection to the GI Genius™ hardware:

- Olympus CV-180 EXERA II;
- Olympus CV-190 EXERA III;
- Olympus CV-1500 EVIS X1;
- Fujifilm VP-4450HD;
- Pentax EPK-i7000 Video Processor;
- Fujifilm VP-7000 ELUXEO.

The GI Genius™ software shall only be used with these VPs.

Medtronic, Medtronic with rising man logo, and Medtronic logo are trademarks of Medtronic. Third-party trademarks (“TM”) belong to their respective owners. The following list includes trademarks or registered trademarks of a Medtronic entity in the United States and/or in other countries:

GI Genius™

Linkverse s.r.l., Via Ostiense 131/L , 00154 Rome, Italy

Distributed by Covidien llc, 15 Hampshire Street, Mansfield, MA 02048 USA.

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
The screenshot shows the 'Status and Settings' interface. On the left, there are three menu items: 'Status' (with a clock icon), 'Volume' (with a speaker icon), and 'Field View' (with a monitor icon). The main area displays 'RX' and 'TX' CPU load gauges at 30% and 40% respectively, and a 'TEMPERATURE' gauge at 45 °C. A QR code and 'CE 0477' are visible at the bottom. Text at the bottom of the screen reads 'Enabled features: detection, characterization v3.0.0 (user manual v3.x)'.

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The screenshot shows the 'Status and Settings' interface with the 'Volume' menu item selected. The main area displays a vertical slider set to 50%. The 'Status' and 'Field View' menu items are visible on the left. Text at the bottom of the screen reads 'Enabled features: detection, characterization v3.0.0 (user manual v3.x)'.

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
The screenshot shows the 'Status and Settings' interface with the 'Volume' menu item selected. The main area displays a vertical slider set to 0%. The 'Status' and 'Field View' menu items are visible on the left. A muted speaker icon is shown next to the 'Volume' menu item. Text at the bottom of the screen reads 'Enabled features: detection, characterization v3.0.0 (user manual v3.x)'.

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The screenshot shows a software interface with a dark background and blue accents. On the left, there is a 'Field View' icon. In the center, there is a QR code with the URL 'linkverse.com/zKb4' below it. On the right, there is a temperature display '60 °C', a 'MD CE 0477' certification mark, and the text 'GI Genius™ Software' followed by the company name 'Linkverse S.r.l.' and address 'Via Ostiense 131L, 00154 Roma, Italy'. At the bottom, there is a status bar with the text 'v3.1.0 - Enabled features: detection, characterization' on the left and '(01)08055060110239(8012)3.1.0' on the right.

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Distribuirano od Covidien LLC, 15 Hampshire Street, Mansfield, MA 02048 SAD.