# 1. Introduction

## 1.1. Safety Symbols

This manual uses the safety symbols below. They denote critical information. Please read them carefully.

|                |  |
|----------------|--|
| <b>DANGER</b>  | Imminently hazardous situation that will result in death or serious injury, if not avoided   |
| <b>WARNING</b> | Potentially hazardous situation that could result in death or serious injury, if not avoided   |
| <b>CAUTION</b> | Potentially hazardous situation that may result in minor or moderate injury, if not avoided.<br>Can alert about unsafe practices or potential equipment damage |
| <b>NOTE</b>    | Useful information   |

## 1.2. Term and abbreviation

The following table summarizes terms and abbreviations used in this document.

| Term | Description            |
|------|------------------------|
| MED  | Main Endoscopy Display |

## 1.3. Document purpose

This document contains essential information on using the USB K-switch in a safe and effective way. Please carefully read this manual and all the equipment manuals involved in the procedure and use the equipment as instructed.

Failure to follow the instructions in this manual may result in damage to and/or malfunction of the equipment. Do not use this device for any other purpose than that for which it has been designed.

This document does not contain any information on how to perform endoscopy procedures, please refer to the related manuals for this information.

Keep this and any related document in a safe and accessible location. If you have questions or comments about this about this manual, please contact helpdesk service (section 7).

## 1.4. Intended use

The USB K-switch is intended to be used by trained clinicians as an accessory to GI Genius™ for the purpose of enabling and disabling the GI Genius™ overlay over Main Endoscopy display (MED).

## 1.5. Device description

The USB K-switch is composed by (Figure 1):

1. Main enclosure;
2. Cable that allows the connection with GI Genius™ hardware;
3. Jack that allows the connection with the Footswitch.

For each press made by the user, the Footswitch sends a switch command to the GI Genius™ hardware through the USB K-switch. The switching command corresponds to the function of the GI Genius™ hardware of switching ON/OFF the software AI of the GI Genius™.

This means that the operator can temporarily halt or restore the GI Genius™ software while in the course of a procedure.

The same functionality is already present on the front panel of the GI Genius™ hardware with a press button, thus the USB K-switch offers an alternative means of using the same functionality.

|               |   |
|---------------|---|
| <b>DANGER</b> | Do not use this USB K-switch for any purpose other than its intended use. |
| <b>DANGER</b> | The USB K-switch does not perform any diagnosis.                          |

The USB K-switch is compliant with IEC/EN 60601-1 and IEC/EN 60601-1-2. This allows the device to be placed in the patient area and to be safely connected to the Footswitch and the GI Genius™ hardware.

## 1.6. Intended users

Users shall have been properly trained on the GI Genius™ device so that they can use and manage it properly.

Users shall be clinicians on lower gastrointestinal endoscopy procedure and the training shall be based on the contents of this user manual.

## 1.7. Symbols and marking

|   |   |
|---|---|
|  Legal Manufacturer   |  Humidity limitation                         |
|  Date of manufacture  |  Temperature limits                          |
|  Do not dispose of the product as normal waste. This is product is subject to Directive 2012/19/EU (WEEE). Follow your local regulations on handling of waste electrical and electronic equipment |  Atmospheric Pressure limitation             |
|  Consult the instruction for use  |  Refer to instruction manual/ booklet (Blue) |
|  The product is a class I medical device in compliance with European Regulation (EU) 2017/745   |  This product is a medical device            |
|  The device is restricted to professional use only  |  Serial Number of the device                 |
|  Keep Dry   |  Catalogue number of the device              |
|  Handle with care   |  Do not use if package is damaged            |
|  Fragile, handle with care  |  Keep away from sunlight                     |
|  Unique Device Identifier   |  Caution                                     |
|   |  Direct Current                              |

Figure 1: USB K-switch general overview

## 2. Information for safe use

### 2.1. USB K-switch compatibility

The USB K-switch is intended to be an accessory to the GI Genius™ hardware and to be connected for such purpose to the Footswitch.

List of compatible devices:

| Part # | Description           |
|--------|-----------------------|
| GG5100 | GI Genius™ System 100 |
| GG5200 | GI Genius™ System 200 |
| GGM100 | GI Genius™ Module 100 |
| GGM200 | GI Genius™ Module 200 |
| FS01   | Footswitch            |

The Footswitch connected to the USB K-switch does operate as a normal SPST (single pole single throw) switch.

In case the Footswitch is not used, it is possible to connect the USB K-switch to an SPST type of switch, either with the CML026R00 cable or directly to the USB K-switch jack port.

Compatible device shall support nominal currents of 10 mA and voltages of 5 V. In case a cable is used for the connection, the maximum resistance shall be 100 Ω.

|               |  |
|---------------|--|
| <b>DANGER</b> | Using incompatible equipment can result in patient injury or equipment damage and makes it impossible to obtain the expected functionality.  |
| <b>DANGER</b> | Connection and installation of this medical electrical equipment with any other equipment shall be responsibility of qualified personnel knowledgeable of electrical safety and electromagnetic interference engineering.<br>Incorrect set up of the end medical system can lead to any of the risks reported in the present manual. |

### 2.2. Installation

|               |   |
|---------------|---|
| <b>DANGER</b> | This device should NEVER be installed or used in areas where the unit could get wet or be exposed to any environmental conditions such as high temperature, humidity, direct sunlight, dust, salt, etc., which could adversely affect the equipment |
| <b>DANGER</b> | This device should NEVER be installed or used in the presence of flammable or explosive gases or chemicals.   |

### 2.3. Repair and modification

|  |   |
|--|---|
| Do not disassemble, modify or attempt to repair the USB K-switch. Any attempt at such operations without the guidance and instructions of technical support will void any warranty of the product. |   |
| <b>DANGER</b>  | Unproper USB K-switch disassembly, modification or repair could result in risks for patients and operators. |

#### 2.3.1. Modification of the setup

|               |   |
|---------------|---|
| <b>DANGER</b> | Setup of the endoscopy column requires technical knowledge on safety aspects, like electrical safety, and operability aspects of the devices, like electromagnetic interferences.<br>Do not attempt to change the setup yourself, but contact the technical department of your Organization or personnel qualified by Linkverse, otherwise injury to patient, operator, or equipment could occur. |
|---------------|---|

### 2.4. Electrical risks

|   |  |
|---|--|
| Following precautions must be strictly observed. Failure to do so may cause danger of electric shock for patient and medical personnel. |  |
| <b>DANGER</b>   | Keep fluids away from the USB K-switch. If fluids are spilled on or into the USB K-switch, immediately stop any operation of the system and contact Support (see section 7).<br>Do not use the USB K-switch when not properly closed or if visually damaged. Do not touch electrical contacts inside any component of the USB K-switch when still connected to GI Genius™. |

### 2.5. Fire or explosion risks

|   |   |
|---|---|
| In order to prevent fire and explosion, do not expose the USB K-switch where: |   |
| <b>DANGER</b>   | there is a high concentration of oxygen.<br>air contains oxidizing agents (e.g. N2O).<br>air contains any flammable gases.<br>there are flammable liquids nearby. |

### 2.6. Interferences and electromagnetic fields

|                |  |
|----------------|--|
| <b>WARNING</b> | When connected to the Footswitch and the GI Genius™ system, the USB K-switch may interfere with other equipment. Please refer to section 1.6 for further information.  |
| <b>CAUTION</b> | Electromagnetic interference may occur to this instrument when it is placed near equipment marked with the following symbol or other portable and mobile RF communications equipment such as mobile phones. If radio interference should occur, mitigation measures may be necessary, such as reorienting or relocating this instrument or shielding the location. Do not use the USB K-switch and its cable closer than 30 cm (12 inches) to the equipment with the following symbol:  |

|                |   |
|----------------|---|
| <b>WARNING</b> | Do not use the USB K-switch where there is a strong electromagnetic field (e.g. Magnetic Resonance, wireless devices, microwaves emitting devices, etc.). |
|----------------|---|

This product is intended for use in the electromagnetic environments specified below. The customer or the user of this product should be sure that it is used in such an environment.

| Electromagnetic emission compliance information and guidance |            |   |
|--|------------|---|
| Emission test  | Compliance | Electromagnetic environment - Guidance  |
| RF emissions<br>CISPR 11                                     | Group 1    | This product uses RF (Radio Frequency) energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.                       |
| RF emissions<br>CISPR 11                                     | Class B    | This product is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |

| Electromagnetic immunity compliance information and guidance   |                             |                             |  |
|--|-----------------------------|-----------------------------|--|
| This product is intended for use in the electromagnetic environments specified below. The customer or the user of this product should assure that it is used in such an environment. |                             |                             |  |
| Immunity test  | IEC 60601-1-2 Test level    | Compliance level            | Electromagnetic environment — Guidance   |
| Electrostatic discharge (ESD)<br>IEC 61000-4-2   | ±8kV contact<br>±15kV air   | ±8kV contact<br>±15kV air   | Floors should be made of wood, concrete, or ceramic tile that hardly produces static. If floors are covered with synthetic material, the relative humidity should be at least 30%. |
| Electrical fast transient/burst<br>IEC 61000-4-4   | ±1kV for input/output lines | ±1kV for input/output lines | Mains power quality should be that of a typical commercial or hospital environment.  |
| Power frequency (50Hz) magnetic field<br>IEC 61000-4-8   | 30A/m                       | 30A/m                       | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.  |
| NOTE UT is the A/C. mains voltage prior to application of the test level.  |                             |                             |  |

| Electromagnetic immunity compliance information and guidance   |                           |                           |  |
|--|---------------------------|---------------------------|--|
| This product is intended for use in the electromagnetic environments specified below. The customer or the user of this product should assure that it is used in such an environment.   |                           |                           |  |
| Immunity test  | IEC 60601-1-2 Test level  | Compliance level          | Electromagnetic environment — Guidance   |
| Conducted RF<br>EN 61000-4-6   | 3Veff<br>150kHz to 80MHz  | 3Veff<br>150kHz to 80MHz  | Portable and mobile RF communications equipment should be used no closer to any part of this product, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.<br>Recommended separation distance<br>$d = 1,2 \times \sqrt{P}$ 150kHz to 80MHz<br>$d = 1,2 \times \sqrt{P}$ 80MHz to 800MHz<br>$d = 2,3 \times \sqrt{P}$ 800MHz to 2,7GHz<br>Where "P" is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and "d" is the recommended separation distance in meters (m). |
| Radiated RF<br>EN 61000-4-3  | 10Veff<br>80MHz to 2,7GHz | 10Veff<br>80MHz to 2,7GHz |  |
| Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:  |                           |                           |  |

## 2.7. Mechanical forces

|                |   |
|----------------|---|
| <b>CAUTION</b> | Do not keep the USB K-switch under the weight of other equipment. They could damage its integrity.                                |
| <b>CAUTION</b> | Once connected to GI Genius™, do not keep the USB K-switch hanging on its cable. It could strain the cable and fall on the floor. |

## 2.8. Connections and disconnections

|                |  |
|----------------|--|
| <b>CAUTION</b> | Do not connect/disconnect the USB K-switch cable while the GI Genius™ is powered ON, otherwise damage to the equipment will occur. |
|----------------|--|

## 2.9. USB K-switch package content

| Part #    | Description              |
|-----------|--------------------------|
| UKS01     | USB K-switch             |
| CML026R00 | USB K-switch Cable       |
| UKS-MN-01 | USB K-switch User Manual |

**WARNING**

Use of accessories other than those specified or provided by the manufacturer of this product could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

### 3. Operating the USB K-switch

#### 3.1. Installation

The USB K-switch is intended to be used connected to the Footswitch and the GI Genius™ hardware.

- Connect the USB K-switch, the Footswitch, and the GI Genius™ hardware.
- Connect the USB cable of USB K-switch with the USB port of GI Genius™ (Figure 2).
- Connect USB K-switch to the Footswitch
- Connect the Footswitch cable with the jack of the USB K-switch (Figure 3).


**Figure 2:** USB K-switch with GI Genius

#### 3.2. Functioning test

Once the devices are correctly connected, follow the steps below in order to verify that the system is correctly working:

- Verify that GI Genius™ is activated and is fully functional: blinking green indicator is shown over the MED;
- Press the Footswitch to disable the GI Genius™ overlay: solid gray strikethrough indicator is shown over the MED;
- Press the Footswitch to disable the GI Genius™ overlay: solid gray strikethrough indicator is shown over the MED;

**Figure 3:** Connection between the Footswitch, the USB K-switch, and the GI Genius hardware

### 4. Technical specifications

|                                   |   |
|-----------------------------------|---|
| Operation environment             | The USB K-switch shall be operating under the following conditions:<br>• $(10 \div 40) ^\circ\text{C} - (50 \div 104) ^\circ\text{F}$<br>• $20\% \div 90\%$ humidity (non-condensing)<br>• $75,3 \text{ kPa} \div 101,3 \text{ kPa}$              |
| Storage and transport environment | The USB K-switch must be stored and transported under the following conditions:<br>• $(-10 \div 40) ^\circ\text{C} - (14 \div 104) ^\circ\text{F}$<br>• $20\% \div 90\%$ humidity (non-condensing)<br>• $75,3 \text{ kPa} \div 101,3 \text{ kPa}$ |
| Power specification               | Supply voltages: 5 VDC / 0,5 A  |

### 5. Maintenance and disposal

Check periodically by visual inspection if the device shows any sign of damage. A damage could lead to the device malfunctioning.

**NOTE**

When disposing of this product, follow all applicable national and local regulations on waste of electrical and electronic equipment.

### 6. Cleaning

Before cleaning, ensure that the GI Genius™ hardware is switched OFF and disconnect the USB K-switch from the GI Genius™ hardware and the Footswitch.

The device does not require regular cleaning, but may be cleaned if dirt/dust has accumulated on visible parts. Wipe clean the device with a soft gauze slightly dampened with water and/or a neutral detergent. Ensure all surfaces are dry before reconnecting to other equipment.

|                |  |
|----------------|--|
| <b>NOTE</b>    | When cleaning the USB K-switch, use only cloth rag dampened with water and a mild cleanser.  |
|                | Do not use cleaning agents which may damage the plastic surfaces, such as detergents, abrasive cleansers or solvent-based cleansers (such as benzene, stain remover).                            |
| <b>WARNING</b> | Do not immerse the USB K-switch into water or other liquids and do not let water droplets fall on the USB port and the female jack port. Water may enter the device and cause electric problems. |

Ensure all surfaces are dry before reconnecting the USB K-switch to the GI Genius™ hardware and the Footswitch. Only clean manually.

### 7. Support

The helpdesk service provided to users of the USB K-switch is available at the following references:

|  |  |
|--|--|
| Email address  | <a href="mailto:sts@cosmopharma.com">sts@cosmopharma.com</a>   |
| <b>NOTE</b>  | In the event that a serious incident has occurred in related to the device use, immediately report the event to the manufacturer, competent authorities, and any other regulators as required. |
|  | Linkverse s.r.l., Via Ostiense 131/L, 00154 Rome, Italy  |

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